

Refrigerated Vaccine Inside at 2°C to 8°C



RECEIVING AND UNPACKING STEPS

- 1 IMMEDIATELY** open the cooler and follow the TagAlert receiver instructions. See below for details.
- 2 Compare** the packing list to the contents of the cooler to confirm they match. See below.
- 3 Inspect** the vaccine and report any tampering (e.g., broken seal). See back for details.
- 4 Place vaccine** in appropriate storage immediately upon receipt.
- 5 Repack** the cooler for refurbishing. See back for details.
- 6 Return** the cooler via UPS. See back for details.



INSTRUCTIONS FOR CONFIRMING PROPER SHIPPING TEMPERATURES

Receiver Instructions

- 1. Stop TagAlert® Indicator**
Press and hold the **START & STOP** button for 5 seconds until Stop icon ● appears.
- 2. Read Indicator Status**
✓ **Accept:** Left arrow ◀ points to Green check mark: Ready for Use.
✗ **Reject:** Right arrow ▶ points to Red X and numbers 1, 2 and/or 3 are visible—a temperature alarm has triggered. **QUARANTINE PRODUCT** and call +1-877-836-7123 or your state/local immunization program right away.

T17003135 Rev. A

TAGALERT INFORMATION

- Remove the TagAlert temperature monitor from the box. The TagAlert is located inside the inner box, next to the vaccine.
- Follow the recipient instructions on the blue card:
 - **Stop TagAlert Indicator.** Press and hold the Start and Stop button until the Stop icon appears.
 - **Read Indicator Status.** Determine if the appropriate temperature was maintained during transit.
 - **TagAlert Monitor Concerns?** Please call **1-877-836-7123**.



Vaccine should be placed in appropriate storage immediately upon receipt.



REVIEW YOUR PACKING SLIP

McKesson Specialty Care Dist. 8745 Forest Hill Irene, Suite 105 Olive Branch, MS 38654		Packing Slip											Page 1 of 2	
		This is not an invoice											Provider PIN: WVA200066 Delivery Number: 230684533 Quality Check Date: 02/11/2025	
Ship-to:		Awarde											Internal use only	
McKesson Specialty Care Dist. 8745 Forest Hill Irene, Suite 105 Olive Branch, MS 38654		McKesson Specialty Care Dist. 8745 Forest Hill Irene, Suite 105 Olive Branch, MS 38654												
NDC	Customer P.O.	Material Description Manufacturer	MFR Lot#	Exp. Date	VFC Doses	317 Doses	State Doses	CHIP Doses	Order Qty	Ship Qty	Unit Price	Extended Price		
19515-0810-52 5017406	394261	FLU; SYR; 10-pack (TIV) PED CDC GLAXO SMITHKLINE	AE2J7	06/16/2025	10	0	0	0	10	10	\$15.22	\$152.20		
Total									10	10		\$152.20		

Verify the **MFR Lot#** and the **Ship Qty** on the packing slip with the contents of the cooler.
The **Ship Qty** is the total number of vaccine cartons.

Refrigerated Vaccine Inside at 2°C to 8°C



WHO TO CALL

Questions?

Need a UPS Pickup of your EcoFlex cooler?

CONTACT

UPS Phone Number: 1-800-377-4877
7:00AM – 9:00PM EST Monday – Friday
7:30AM – 6:00PM EST Saturday

UPS Email: Healthcaresupport@UPS.com
7:00AM – 9:00PM EST Monday – Friday
7:30AM – 6:00PM EST Saturday


WEB: UPS.com/healthcare

Provide the McKesson Account Number: 20V8R9



RETURN YOUR COOLER

PLEASE
RETURN COOLER
DO NOT
DISPOSE



Step 1: Open shipper and remove the top PCM gel pack(s).



Step 2: Remove the product(s) from the inner corrugated box.



Step 3: Ensure the inner corrugated box is back in the EcoFlex shipping system.



Step 4: Place the top PCM gel pack(s) back on top.



Step 5: Fold in the original outer flaps. Fold in the remaining flaps with return.



Step 6: Schedule a pickup with UPS. Contact UPS either by phone or by email.



The return shipping label can be found on the inside flap of your EcoFlex box. For more information, contact info@coldchaintech.com

McKesson